

Specifications

Product name	Konan Specular Microscope XX
Model no.	CellChek 20-1 CellChek 20-1 PLUS
Type	Class I, Type B
Photographic field	0.25 x 0.55mm
Photographic capability	Automatic
Photographic location	20-1: Center, 6 points in the periphery (12, 2, 4, 6, 8, 10 o' clock) 20-1 PLUS: Center, 6 points in the periphery (12, 2, 4, 6, 8, 10 o' clock), 6 points in the outer periphery (2, 3, 4, 8, 9, 10 o' clock)
Central cornea thickness	Measurable
Illumination	LED
Analytical method	Auto Center Method / Auto F Center Method / Auto Trace Method / Center Method / F Center Method / Trace Method
Analysis Data	Cell density (mm ²), Coefficient of Variation, Standard Deviation, Percentage of hexagonal cells, Numbers of calculated cells, Average cell area (μm ²), Maximum cell area (μm ²),
Integrated monitor	10.1" Capacitive touch panel with 180 degrees flexible horizontal & vertical movement capability.
External interface	EMR connection / DICOM capability
External ports	USB3.0 x 4, RJ45
Power	100-240 VAC, 50/60 Hz

N 0123



The All New, *Blazing Fast, Fully-Automated* Specular Microscope



Interested in learning more?

- +1 (949) 576-2200
- Live Chat @ KonanMedical.com
- Sales@KonanMedical.com
- KonanMedical.com

1. Based on benchmark data..



Specular Microscopy

THE GOLD STANDARD REDEFINED.

Faster, Easier Endothelial Imaging

“Process bilateral exams with one touch in under 40 seconds.”

Introducing the all new CellChek® 20

CellChek 20 is Konan Medical’s new, non-contact specular microscope that can capture and analyze bilateral exams with one touch in *under 40 seconds*¹.

The new, easy-to-use “Simple Mode” enables one-touch, fully-automated endothelial image capture with analysis, reporting, and exporting of data.

The completely reimagined CellChek 20 is noticeably improved and simplified; *remove the device from the box, connect to power, turn it on, and start imaging.*

Konan’s gold standard *Center Method™* and *Flex Center Method™* are now *fully-automated*. Automated Center Method displays cell center counted densities and morphometric indices of both eyes *immediately after* images are captured – *all in under 40 seconds*.



The flexible 10.6” wide touch screen can be turned and tilted 180°, providing better ergonomics, space savings, and flexibility with technician use from any of the four sides.

CAPTURE, ANALYZE, & OUTPUT EXAMS FOR BOTH EYES, WITH JUST ONE TOUCH.

Key New Features

With CellChek 20 you can capture, analyze and print/export a bilateral exam with just one touch. The new, easy-to-use “Simple Mode” makes running a bilateral exam faster and easier than ever before. CellChek 20 reduces testing time significantly, allowing you to test patients quickly and efficiently.

One-Touch Imaging

Fast, one-touch imaging of *both eyes* including positioning, alignment, imaging capture, analysis and print-out/exporting data.

Optimized Touch Functions

Easy touch-friendly interface makes capturing, analysis, and reporting fast and efficient. Also features pinch-to-zoom image adjustment.

Automated Trends Analysis

Konan exclusive image location enables relevant visit-to-visit trend analysis to track patient progression.

25% Higher Resolution

Improved resolution allows better visualization of single cells and cell walls. Superb cell border tracing and gutta exclusion.

Auto-Capture/Retry

This improved technology allows for successful imaging of difficult or challenging cases.

> 30% Larger Image Area


The field of view has been increased to 0.25mm x 0.55mm, enabling better corneal visualization and more cells to be analyzed.

Fully-Automated Analysis

CellChek 20 offers *fully-automated*: Center Method; Flex-Center Method; and Auto-Trace analysis. These analysis methods provide one-touch, high precision results *FAST*.

These fully-automated analysis methods are offered alongside the traditional manual Center and Flex-Center analysis methods, offering a total of 6 analysis options.

Konan’s Center Method is mentioned in FDA panel minutes as being the “gold standard”, and is used by virtually every professional reading center.



Auto-Center Method
With the Auto-Center Method, the center of each individual, visible cell is automatically marked. Guttata and other dark regions are automatically excluded.

Auto Flex-Center Method
The Auto Flex-Center Method is an additional automated tool for analyzing images of advanced stage diseased corneas in which only a very few cells are visible.

